


The SVR System


Create an Elliptical Ventricle Every Time

The SVR® procedure is a comprehensive surgical therapy for patients suffering from ischemic heart failure. Studies have shown that surgically restoring a dilated left ventricle to a more normal elliptical anatomy using the SVR system results in significantly improved ejection fraction and less risk of post-operative mitral regurgitation.

Elliptical Ventricle Post SVR
Optimal surgical outcome



Spherical Ventricle Post SVR
Undesirable surgical outcome



An SVR technique should result in a normal elliptical ventricle every time.

Key Benefits:

- Standardizes the SVR procedure
- Achieves optimal ventricular function and volume
- Defines the new apex
- Reorients the mitral apparatus and papillary muscles properly

The SVR System

Standardizing Surgical Ventricular Restoration (SVR) for consistent and elliptical results every time.

The SVR System includes:

The Mannequin™


- Defines the new apex
- Acts as a template to create an elliptical ventricle with optimal volume every time
- Offered in various sizes based on the patient's BSA

Patch Sizer

- Reduces potential for incorrectly trimming the patch
- Decreases the procedure time

Also included are conveniently sized syringes and a marking pen.

Ordering Information

Product*	Code	Description	Units
	SVC-0080S	SVR System with 80cc Mannequin, without patch	1
	SVC-0085S	SVR System with 85cc Mannequin, without patch	1
	SVC-0090S	SVR System with 90cc Mannequin, without patch	1
	SVC-0100S	SVR System with 100cc Mannequin, without patch	1
	SVC-0110S	SVR System with 110cc Mannequin, without patch	1
	SVC-0120S	SVR System with 120cc Mannequin, without patch	1
	SVC-0130S	SVR System with 130cc Mannequin, without patch	1
	SVC-0140S	SVR System with 140cc Mannequin, without patch	1
	SVC-0150S	SVR System with 150cc Mannequin, without patch	1

* Patents US6681773 and additional patents pending

* Dacron Hemashield Patch sold separately

To Place Orders

Chase Medical Customer Service

Ph: 972-783-7005 or 1-800-787-0378

Fx: 972-235-3446 or 1-844-272-6624

Email: cust_svc@chasemedical.com

Chase Medical, 885 East Collins Boulevard, Suite 110, Richardson, Texas 75081 U.S.A. / www.chasemedical.com